

Anti-Human/Mouse/Rat EPHA3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03765

3th Edition

Description: Receptor for members of the ephrin-A family. Binds to ephrin-A2, -A3, -A4 and -A5. Could play a role in lymphoid function.

Antigen: Synthesized peptide derived from human EphA3 around the non-phosphorylation site of Y602.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed. Highest level in placenta.