

Anti-Human/Mouse/Rat Ephrin A5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03792

3th Edition

Description: May function actively to stimulate axon fasciculation. Induces compartmentalized signaling within a caveolae-like membrane microdomain when bound to the extracellular domain of its cognate receptor. This signaling event requires the activity of the Fyn tyrosine kinase.

Antigen: Synthesized peptide derived from the Internal region of human Ephrin-A5

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: