

**Anti-Human/Mouse/Rat ERK 3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA03824**

---

3th Edition

**Description:** Phosphorylates microtubule-associated protein 2 (MAP2). May promote entry in the cell cycle.

**Antigen:** Synthesized peptide derived from the Internal region of human ERK 3

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Highest expression in the skeletal muscle, followed by the brain. Also found in heart, placenta, lung, liver, pancreas, kidney and skin fibroblasts.