

Anti-Human/Mouse ERK 8 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03841

3th Edition

Description:In vitro, phosphorylates MBP.

Antigen:Synthesized peptide derived from human ERK 8 around the non-phosphorylation site of Thr175/Tyr177.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Widely expressed with a maximal expression in lung and kidney.