

Anti-Human/Mouse/Rat TRIM72 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03880

3th Edition

Description:TRIM72 (Tripartite Motif Containing 72) is a Protein Coding gene. Among its related pathways are Cardiac conduction and Smooth Muscle Contraction. GO annotations related to this gene include phosphatidylserine binding. An important paralog of this gene is TRIM50.

Antigen:Recombinant Protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: