

Anti-Human/Mouse/Rat FER Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03912

3th Edition

Description: Tyrosine kinase of the non-receptor type. Probably performs an important function, perhaps in regulatory processes such as cell cycle control.

Antigen: Synthesized peptide derived from the C-terminal region of human Fer

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in a variety of lymphoid cell lines.