

Anti-Human/Mouse FGRL1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03922

3th Edition

Description:Has a negative effect on cell proliferation.

Antigen:Synthesized peptide derived from the Internal region of human FGFR-5

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed preferentially in cartilaginous tissues and pancreas. Highly expressed in the liver, kidney, heart, brain and skeletal muscle. Weakly expressed in the lung, small intestine and spleen.