

Anti-Human/Mouse FARP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03931

3th Edition

Description:FARP2 (FERM, ARH/RhoGEF And Pleckstrin Domain Protein 2) is a Protein Coding gene. Diseases associated with FARP2 include Pseudohypoparathyroidism Ia and Aarskog-Scott Syndrome. Among its related pathways are Development Endothelin-1/EDNRA signaling and Focal Adhesion. GO annotations related to this gene include Rho guanyl-nucleotide exchange factor activity and Rac guanyl-nucleotide exchange factor activity. An important paralog of this gene is FARP1.

Antigen:Synthesized peptide derived from the Internal region of human FIR

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: