

Anti-Human/Mouse/Rat FST Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03956

3th Edition

Description: Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

Antigen: Synthesized peptide derived from the Internal region of human Follistatin

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Isoform 1 is the predominant isoform in serum but is undetectable in follicular fluid.