

Anti-Human/Mouse/Rat FGFR1OP Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03965

3th Edition

Description: Required for anchoring microtubules to the centrosomes.

Antigen: Synthesized peptide derived from the C-terminal region of human FOP

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous. Highly expressed in heart, liver, muscle, kidney, intestine, colon, adrenal gland, prostate, testis, and pancreas.