

Anti-Human/Mouse/Rat IL22 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03991

3th Edition

Description:IL22 (Interleukin 22) is a Protein Coding gene. Diseases associated with IL22 include Acute Generalized Exanthematous Pustulosis and Chronic Mucocutaneous Candidiasis. Among its related pathways are PEDF Induced Signaling and Common Cytokine Receptor Gamma-Chain Family Signaling Pathways. GO annotations related to this gene include cytokine activity and interleukin-22 receptor binding.

Antigen:Synthesized peptide derived from the C-terminal region of human IL-22.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: