

Anti-Human/Mouse DROSHA Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA03997

3th Edition

Description: This gene encodes a ribonuclease (RNase) III double-stranded RNA-specific ribonuclease and subunit of the microprocessor protein complex, which catalyzes the initial processing step of microRNA (miRNA) synthesis. The encoded protein cleaves the stem loop structure from the primary microRNA (pri-miRNA) in the nucleus, yielding the precursor miRNA (pre-miRNA), which is then exported to the cytoplasm for further processing. In a human cell line lacking a functional copy of this gene, canonical miRNA synthesis is reduced. Somatic mutations in this gene have been observed in human patients with kidney cancer. DROSHA (Drosha Ribonuclease III) is a Protein Coding gene. Diseases associated with DROSHA include Combat Disorder and Murray Valley Encephalitis. Among its related pathways are mRNA Capping and Gene Expression. GO annotations related to this gene include poly(A) RNA binding and RNA binding.

Antigen: Synthesized peptide derived from the Internal region of human RNase III Drosha.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous.