

Anti-Human/Mouse/Rat APOOL Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04045

3th Edition

Description:APOOL (Apolipoprotein O Like) is a Protein Coding gene. Among its related pathways are Response to elevated platelet cytosolic Ca²⁺. An important paralog of this gene is APOO.This gene encodes a protein which contains an apolipoprotein O superfamily domain. This domain is found on proteins in circulating lipoprotein complexes.

Antigen:Synthesized peptide derived from the C-terminal region of human ApoOL

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: