

**Anti-Human/Mouse/Rat FOXE1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04064**

---

3th Edition

**Description:** Probable transcription factor. Could be involved in thyroid gland organogenesis.

**Antigen:** Synthesized peptide derived from the N-terminal region of human FoxE1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Detected in adult brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, heart, colon, small intestine testis and thymus. Expression was strongest in heart and pancreas.