

Anti-Human/Mouse FRS2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04095

3th Edition

Description: Adapter protein that links FGR and NGF receptors to downstream signaling pathways. Involved in the activation of MAP kinases. Modulates signaling via SHC1 by competing for a common binding site on NTRK1.

Antigen: Synthesized peptide derived from human FRS2 around the non-phosphorylation site of Y196.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in heart, brain, spleen, lung, liver, skeletal muscle, kidney and testis.