

Anti-Human/Monkey FRS3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04097

3th Edition

Description: SNT2 is a substrate for the fibroblast growth factor receptor. It functions in linking FGF receptor stimulation to activators of Ras.

Antigen: Synthesized peptide derived from the Internal region of human FRS3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: