



**Anti-Human/Mouse/Rat NOLC1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04107**

---

3th Edition

**Description:**NOLC1 (Nucleolar And Coiled-Body Phosphoprotein 1) is a Protein Coding gene. Diseases associated with NOLC1 include Treacher Collins Syndrome 1. GO annotations related to this gene include poly(A) RNA binding and GTP binding.

**Antigen:**Synthesized peptide derived from the C-terminal region of human Nopp140

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**