

Anti-Human/Mouse/Rat G3BP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04108

3th Edition

Description:G3BP2 (G3BP Stress Granule Assembly Factor 2) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is G3BP1.

Antigen:Synthesized peptide derived from the C-terminal region of human G3BP2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: