

Anti-Human/Mouse/Rat GALK1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04128

3th Edition

Description:Galactokinase is a major enzyme for the metabolism of galactose and its deficiency causes congenital cataracts during infancy and presenile cataracts in the adult population.

Antigen:Synthesized peptide derived from the N-terminal region of human GALK1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: