

Anti-Human GSDMC Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04135

3th Edition

Description:GSDMC (Gasdermin C) is a Protein Coding gene. An important paralog of this gene is GSDMA.

Antigen:Synthesized peptide derived from the Internal region of human Gasdermin C

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed mainly in trachea and spleen (PubMed:11223543). In the esophagus, expressed in differentiating cells and probably in differentiated cells. Also detected in gastric epithelium (PubMed:19051310).