

Anti-Human/Mouse/Rat EPG5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04141

3th Edition

Description: This gene encodes a large coiled coil domain-containing protein that functions in autophagy during starvation conditions. Mutations in this gene cause Vici syndrome. EPG5 (Ectopic P-Granules Autophagy Protein 5 Homolog) is a Protein Coding gene. Diseases associated with EPG5 include Vici Syndrome and Spastic Paraplegia 49, Autosomal Recessive.

Antigen: Recombinant Protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: