

Anti-Human GBP3 Polyclonal Antibody
Polyclonal Antibody

Cat.NO.: PA04148

3th Edition

Description: This gene encodes a member of the guanylate-binding protein (GBP) family. GBPs specifically bind guanine nucleotides (GMP, GDP, and GTP) and contain two of the three consensus motifs found in typical GTP-binding proteins. The encoded protein interacts with a member of the germinal center kinase family. Multiple transcript variants encoding different isoforms have been found for this gene. GBP3 (Guanylate Binding Protein 3) is a Protein Coding gene. Among its related pathways are Interferon gamma signaling and Innate Immune System. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is GBP1.

Antigen: Synthesized peptide derived from the C-terminal region of human GBP3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: