

Anti-Human GBP4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04150

3th Edition

Description: Guanylate-binding proteins, such as GBP4, are induced by interferon and hydrolyze GTP to both GDP and GMP. GBP4 (Guanylate Binding Protein 4) is a Protein Coding gene. Among its related pathways are Interferon gamma signaling and Innate Immune System. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is GBP7.

Antigen: Synthesized peptide derived from the C-terminal region of human GBP4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: