

Anti-Human/Mouse/Rat TLR1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04156

3th Edition

Description: Participates in the innate immune response to microbial agents. Specifically recognizes diacylated and triacylated lipopeptides. Cooperates with TLR2 to mediate the innate immune response to bacterial lipoproteins or lipopeptides (PubMed:21078852). Forms

Antigen: Recombinant Protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous. Highly expressed in spleen, ovary, peripheral blood leukocytes, thymus and small intestine.