

**Anti-Human/Mouse/Rat TLR1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04158**

---

3th Edition

**Description:** Participates in the innate immune response to microbial agents. Specifically recognizes diacylated and triacylated lipopeptides. Cooperates with TLR2 to mediate the innate immune response to bacterial lipoproteins or lipopeptides (PubMed:21078852). Forms

**Antigen:** Recombinant Protein

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous. Highly expressed in spleen, ovary, peripheral blood leukocytes, thymus and small intestine.