

Anti-Human/Mouse GCP3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04159

3th Edition

Description:TUBGCP3 (Tubulin Gamma Complex Associated Protein 3) is a Protein Coding gene. Diseases associated with TUBGCP3 include Cardiomyopathy, Dilated, 1P. Among its related pathways are Regulation of PLK1 Activity at G2/M Transition and Chks in Checkpoint Regulation. GO annotations related to this gene include structural molecule activity and gamma-tubulin binding. An important paralog of this gene is TUBGCP2.

Antigen:Synthesized peptide derived from the Internal region of human GCP3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitously expressed.