

**Anti-Human/Mouse/Rat GCP5 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04163**

---

3th Edition

**Description:**TUBGCP5 (Tubulin Gamma Complex Associated Protein 5) is a Protein Coding gene. Among its related pathways are Regulation of PLK1 Activity at G2/M Transition and Chks in Checkpoint Regulation. GO annotations related to this gene include microtubule binding. An important paralog of this gene is TUBGCP6.

**Antigen:**Synthesized peptide derived from the Internal region of human GCP5

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Widely expressed, with highest levels in heart and skeletal muscle and moderate levels in brain.