

**Anti-Human/Mouse/Rat STAT2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04183**

---

3th Edition

**Description:** Signal transducer and activator of transcription that mediates signaling by type I IFNs (IFN-alpha and IFN-beta). Following type I IFN binding to cell surface receptors, Jak kinases (TYK2 and JAK1) are activated, leading to tyrosine phosphorylation of STA

**Antigen:** Synthesized peptide derived from human Stat2 around the non-phosphorylation site of Y631.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**