

Anti-Human/Mouse/Rat GIT1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04194

3th Edition

Description: GIT1 (GIT ArfGAP 1) is a Protein Coding gene. Diseases associated with GIT1 include Chromosome 17Q11.2 Deletion Syndrome, 1.4-Mb and Chromosome 17Q11.2 Deletion Syndrome. Among its related pathways are Arf6 signaling events and Adhesion. GO annotations related to this gene include GTPase activator activity. An important paralog of this gene is GIT2.

Antigen: Synthesized peptide derived from the C-terminal region of human GIT1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: