

Anti-Human/Mouse/Rat Glucagon Receptor Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04210

3th Edition

Description:The protein encoded by this gene is a glucagon receptor that is important in controlling blood glucose levels. Defects in this gene are a cause of non-insulin-dependent diabetes mellitus (NIDDM).

Antigen:Synthesized peptide derived from the Internal region of human Glucagon Receptor

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: