

Anti-Human STAT5A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04220

3th Edition

Description: Carries out a dual function: signal transduction and activation of transcription. Binds to the GAS element and activates PRL-induced transcription.

Antigen: Synthesized peptide derived from human Stat5a around the non-phosphorylation site of S780.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: