

Anti-Human GPR113 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04247

3th Edition

Description:GPR113 belongs to the G-protein coupled receptor 2 family, LN-TM7 subfamily. It contains one GPS domain. GPR113 is an orphan receptor.

Antigen:Synthesized peptide derived from the N-terminal region of human GPR113

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: