

Anti-Human/Mouse/Rat/Pig IGFN1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04251

3th Edition

Description:IGFN1 (Immunoglobulin-Like And Fibronectin Type III Domain Containing 1) is a Protein Coding gene. An important paralog of this gene is IGSF22.

Antigen:Recombinant Protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in skeletal muscle.