

Anti-Human/Mouse/Rat GPR139 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04265

3th Edition

Description:Orphan receptor. Seems to act through a G(q/11)-mediated pathway.

Antigen:Synthesized peptide derived from the Internal region of human GPR139

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed almost exclusively in the brain. Detected at very low levels in the peripheral tissues.