

**Anti-Human/Mouse/Rat GPR173 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04281**

---

3th Edition

**Description:**Orphan receptor.

**Antigen:**Synthesized peptide derived from the Internal region of human GPR173

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed at high levels in brain and ovary. Lower levels in small intestine. In brain regions, detected in all regions tested. Highest levels in the cerebellum and cerebral cortex.