

Anti-Human/Mouse/Rat GPR175 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04283

3th Edition

Description:TPRA1 (Transmembrane Protein Adipocyte Associated 1) is a Protein Coding gene. GO annotations related to this gene include G-protein coupled receptor activity.

Antigen:Synthesized peptide derived from the Internal region of human GPR175

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous, with higher levels in heart, placenta and kidney.