

**Anti-Human GPR92 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04315**

---

3th Edition

**Description:** Receptor for lysophosphatidic acid (LPA), a mediator of diverse cellular activities.

**Antigen:** Synthesized peptide derived from the Internal region of human GPR92

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Not expressed in frontal cortex, basal forebrain, caudate putamen, thalamus, or hippocampus.