

Anti-Human/Mouse/Rat KCNA3 Polyclonal Antibody**Polyclonal Antibody****Cat.NO.: PA04320**

3th Edition

Description: Potassium channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. Four sequence-related potassium channel genes - shaker, shaw, shab, and shal - have been identified in *Drosophila*, and each has been shown to have human homolog(s). This gene encodes a member of the potassium channel, voltage-gated, shaker-related subfamily. This member contains six membrane-spanning domains with a shaker-type repeat in the fourth segment. It belongs to the delayed rectifier class, members of which allow nerve cells to efficiently repolarize following an action potential. It plays an essential role in T-cell proliferation and activation. This gene appears to be intronless and it is clustered together with KCNA2 and KCNA10 genes on chromosome 1. KCNA3 (Potassium Voltage-Gated Channel Subfamily A Member 3) is a Protein Coding gene. Diseases associated with KCNA3 include Lymph Node Cancer. Among its related pathways are Potassium Channels and Transmission across Chemical Synapses. GO annotations related to this gene include ion channel activity and voltage-gated ion channel activity. An important paralog of this gene is KCNA2.

Antigen: Synthetic Peptide**Form:****How to use:** 1.0 ml distilled water will be added to the product**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C**Dilution:** PBS (pH7.4) containing 1% BSA**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.**Specificity:**