



**Anti-Human GPRIN2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04335**

---

3th Edition

**Description:**GPRIN2 (G Protein Regulated Inducer Of Neurite Outgrowth 2) is a Protein Coding gene. An important paralog of this gene is GPRIN1.

**Antigen:**Synthesized peptide derived from the N-terminal region of human GRIN2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**