

Anti-Human/Mouse/Rat GRK5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04347

3th Edition

Description: Specifically phosphorylates the activated forms of G protein-coupled receptors.

Antigen: Synthesized peptide derived from the Internal region of human GRK 5

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highest levels in heart, placenta, lung > skeletal muscle > brain, liver, pancreas > kidney.