



**Anti-Human TCEAL6 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04348**

---

3th Edition

**Description:** TCEAL6 (Transcription Elongation Factor A Like 6) is a Protein Coding gene. An important paralog of this gene is TCEAL3.

**Antigen:** Synthesized peptide derived from the Internal region of human TCEAL6

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**