

Anti-Human GRK7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04350

3th Edition

Description: Phosphorylates cone opsins thereby initiating their deactivation.

Antigen: Synthesized peptide derived from the Internal region of human GRK 7

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Retinal cones, outer and inner segments.