

Anti-Human GRTP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04352

3th Edition

Description:GRTP1 (Growth Hormone Regulated TBC Protein 1) is a Protein Coding gene. GO annotations related to this gene include GTPase activator activity. An important paralog of this gene is TBC1D10C.

Antigen:Synthesized peptide derived from the Internal region of human GRTP1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: