

**Anti-Human ANKRD20A1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04355**

---

3th Edition

**Description:**ANKRD20A1 (Ankyrin Repeat Domain 20 Family Member A1) is a Protein Coding gene. Diseases associated with ANKRD20A1 include Behcet Syndrome. An important paralog of this gene is ANKRD20A2.

**Antigen:**Synthesized peptide derived from the Internal region of human A20.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**