

**Anti-Human Histone 1.0 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04391**

---

3th Edition

**Description:** Histones H1 are necessary for the condensation of nucleosome chains into higher-order structures. The H1F0 histones are found in cells that are in terminal stages of differentiation or that have low rates of cell division.

**Antigen:** Synthesized peptide derived from the Internal region of human Histone 1.0

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**