

**Anti-Human/Mouse/Rat/Monkey Histone deacetylase 10 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04392**

---

3th Edition

**Description:**Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes.

**Antigen:**Synthesized peptide derived from the N-terminal region of human Histone deacetylase 10

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitous. High expression in liver, spleen, pancreas and kidney.