

Anti-Human/Mouse HNRNPA2B1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04442

3th Edition

Description:Involved with pre-mRNA processing. Forms complexes (ribonucleosomes) with at least 20 other different hnRNP and heterogeneous nuclear RNA in the nucleus.

Antigen:Synthesized peptide derived from the N-terminal region of human hnRNP A2/B1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: