

**Anti-Human HOXA6 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04455**

---

3th Edition

**Description:**In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which Human, Mouse regulate gene expression, morphogenesis, and differentiation.

**Antigen:**Synthesized peptide derived from the Internal region of human HoxA6

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**