

Anti-Human/Mouse/Rat SSTR2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA04459

3th Edition

Description: Somatostatin acts at many sites to inhibit the release of many hormones and other secretory proteins. The biologic effects of somatostatin are probably mediated by a family of G protein-coupled receptors that are expressed in a tissue-specific manner. SSTR2 is a member of the superfamily of receptors having seven transmembrane segments and is expressed in highest levels in cerebrum and kidney.

Antigen: Synthesized peptide derived from the Internal region of human SSTR2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Cerebrum and kidney. In lesser amounts in jejunum, colon and liver.