

**Anti-Human/Mouse/Rat MUM1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA04480**

---

3th Edition

**Description:**MUM1,also named as EXPAND1,is an architectural component of the chromatin,which in response to DNA damage serves as an accessory factor to promote cell survival. It is a 53BP1-associated DNA damage-responsive factor. MUM1/EXPAND1 is different to another protein IRF4/MUM1.

**Antigen:**Synthesized peptide derived from MUM1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**